

AUR-400 URINE ANALYZER

AUR-400 Urine Analyzer is a device designed for early screening and detecting kidney diseases, diabetes, urinary tract infection and liver diseases.

Advantages:

- ◆ Comprehensive Range of Test Menu
- ◆ Testing Up to 14 Parameters in One Go
- ◆ Language Customizable
- ◆ Built-in Thermal Printer for Instant Printing
- ◆ Analyte Combination Customizable
- ◆ External WIFI Connection Function Optional
- ◆ Touch Screen Operation with Big Screen for Easy Operation

**TECHNICAL SPECIFICATIONS**

Model	AUR-400
Principle	Photoelectric Colorimetry
Testing Speed	600 tests/hour
Data Memory	9999 results
Strip Incubation Time	1 minute
Light Wavelength	460nm, 550nm and 620nm
Display	7" LCD touch screen
Printer	Built-in thermal printer
Connectivity	LAN port, USB 2 port, Type-C port, COM port, External WIFI (Optional)
Compatibility	LIS
Working Condition	10-30°C, RH≤85%, 86-106kPa
Operating Voltage	AC100-240V, 50/60Hz
Weight	6.5 kg
Dimension	380*330*215mm