

TD3 TABLE TOP LOW SPEED CENTRIFUGE

- ◆ Microprocessor control, DC brushless motor for stable and low noise operation.
- ◆ Parameters can be set freely during a run
- ◆ Countdown time showed by seconds if the centrifugation time is less than 1min.
- ◆ Self-diagnosis system can automatically detect the faults of overspeed and no speed signal, door lid interlock to improve the safety performance.
- ◆ Compact with great performance.
- ◆ Special damping structure to reduce vibration and noise.

**TECHNICAL PARAMETER**

Model	TD3
Max. Speed	4000rpm
Max. RCF	1980xg
Max. Capacity	8 x 15ml
Speed Accuracy	±30rpm
Time Setting Range	1min to 99min
Noise	<62dB(A)
Power Supply	AC220V±22V 50/60Hz 2A
Total Power	100W
Dimension (W x D x H)	265x305x250mm
Package Dimension (W x D x H)	322x452x268mm
Net Weight	7.5kg



Capacity: 8x15ml
Max. Speed: 4000 rpm
Max. RCF: 1980xg