

L500 TABLETOP LOW SPEED CENTRIFUGE

L500 Tabletop Low Speed Centrifuge is designed for the clinical research with our years of experience. Stainless rotors with variety of capacity for convenient operation.

Features:

- ◆ Microprocessor control and DC brushless motor.
- ◆ Touch panel and LCD display.
- ◆ RCF value can be calculated automatically.
- ◆ Special damping structure to reduce vibration.
- ◆ Electrical door interlock, with which centrifuge cannot work if the door is open and the door cannot be opened when it is working.
- ◆ Variety of brackets for users to convenient choice.

**TECHNICAL PARAMETERS**

Model	L500
Max. Speed	5000rpm
Max. RCF	4030xg
Max. Capacity	16x15ml
Speed Accuracy	±30rpm
Time Setting Range	1min to 99min
Noise	≤65dB(A)
Power Supply	AC220V±22V 50/60Hz 5A
Total Power	250W
Dimension (W x D x H)	370x432x250mm
Package Dimension(W x D x H)	435x580x312mm
Net Weight	15kg



Capacity: 24x10ml
Max. Speed: 5000rpm
Max. RCF: 4030xg



Capacity: 4x50ml, 4x4x15ml, 4x6x5ml
Max. Speed: 4200 rpm
Max. RCF: 2760xg