


LT53 TABLETOP LOW SPEED CENTRIFUGE

LT53 Tabletop Low Speed Centrifuge is widely used in clinics, genetic biology and cytology with its great performance and high quality. Various rotors, brackets and adaptors are available for choice.

Features:

- ◆ Microprocessor control, DC brushless motor.
- ◆ Electronic door interlock for safe operation.
- ◆ RCF value can be calculated automatically.
- ◆ Touch panel and LCD display for easy operation.
- ◆ Various rotors, brackets and adaptors for choice.
- ◆ It is programmable and 25 programs can be stored.



TECHNICAL PARAMETER

Model	LT53
Max. Speed	5300rpm
Max. RCF	4113xg
Max. Capacity	4x250ml
Speed Accuracy	±20rpm
Time Setting Range	1min to 99min
Noise	≤65dB(A)
Power Supply	AC220V±22V 50/60Hz
Total Power	400W
Dimension (W x D x H)	436x560x380mm
Package Dimension (W x D x H)	547x727x432mm
Net Weight	42kg



Capacity: 4x100ml



Capacity: 4x2x50ml, 4x2x100ml
4x8x15ml, 4x16x5ml/7ml
4x12x5ml/7ml



Capacity: 4x250ml



Capacity: 2x2x96 well



Capacity: 48x10ml