

URINANALYSIS REAGENT STRIPS

The Urinalysis Reagent Strips (Urine) are firm plastic strips onto which several separate reagent areas are affixed. The test is for the qualitative and semi-quantitative detection of one or more of the following analytes in urine: Ascorbic acid, Creatinine, Microalbumin, Glucose, Bilirubin, Ketone Acetoacetic acid, Specific Gravity, Blood, pH, Protein, Urobilinogen, Nitrite and Leukocytes, Calcium. The Urinalysis Reagent Strips (Urine) are for single use in professional near-patient (point-of-care) and centralized laboratory locations.

Specifications:

- ◆ For rapid detection of multiple analytes in human urine.
- ◆ For professional in vitro diagnostic use only.
- ◆ Do not use after the expiration date.
- ◆ Do not re-use.
- ◆ Suitable for both reading by eye and device (urine analyzer).

